

# TAPPED DENSITY TESTER

303-D



## Product description

- Equipment designed for powder, granules and similar product's density measurement and tests, in conformity with USP Method 1 and Method 2
- Extern structure: acoustic cabinet to prevent loud noises in the lab area and glass window.

## Technical Characteristics

- LCD display control panel with Membrane function keypad for tap/min and test-run time selection
- Tap Count: 250 to 300 RPM selectable on the equipment
- Drop height: 14.0mm to 3.0mm, selectable on the equipment.

## Technical Specifications

Timer	1 – 600 m
Tap count	1 - 9999
Electrical specs	110V or 220V 50Hz/60Hz
Protection grade	IP 20

## INCLUDED ACCESSORIES AND COMPONENTS



- Operational manual
- 250 ml and 100 ml glass cylinders

## OPTIONAL



- IQ and OQ qualification protocols

## Taped density tester 303D information

WIDTH mm	DEPTH mm	HEIGHT mm	ELECTRICAL SPECS W	WEIGHT Kg
560	510	710	100	37